

Product Name : Transport Incubator	Product Code : NLABMEDICAL-M2740007
Description : Transport Incubator	
Technical Specification : Be it inter or intra hospital transfer, proper transportation is crucial from the time infant is picked up till reaching the destination. Microprocessor controlled system LED display Double wall made of high optical quality acrylic Precise distribution of heat over the whole hood Removable parts easy to clean Internal Noise Level <60 dB Memory foam mattress for newborn comfort Four Oval Ports Easy access panel in front and back Diffuser for Pipe Passage, 2 front & 2 back Patient skin temperature probe Passive Humidity System Skin mode and Air control mode Air Temperature monitoring Patient temperature monitoring	

Oxygen concentration

Comprehensive alarm system

Heating failure alarm

Skin probe failure

Air probe failure

Air circulation failure

Power supply failure

System failure

Low battery

Air micro-filter

Castors with brakes

IV Pole

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003