

Product Name :
Transport Ventilator

Product Code :
NLABMEDICAL-M3390001



Description :

Transport Ventilator

Technical Specification :

Dimensions W 50 X D 23 X H 23 cm (basic dimensions)

Weight Approximately 19 Kg (42 lb)

Insp. flow Constant or accelerating.

Pattern decelerating in Pressure

Controlled Ventilation

Compressible 0.075 l (internal)

Working pressure Up to 120 cm H₂O

Preset inspiratory 0.5-40 l/in. With sufficient working minute volume pressure

Respiratory rate 5-120 breaths/min

Inspiratory time 20, 25, 33, 50, 67 or 80 % of the respiratory cycle

Pause time 0, 5, 10, 20 or 30 % of the respiratory cycle

PEEP 0-50 cm H₂O

Upper Pressure 15-120 cm H₂O limit

Trig Sensitivity 20-0 cm H₂O (below PEEP)

SIMV-rate HIGH RATE (4-40 SIMV breaths/min)

Low rate (0.4 -40SIMV breaths/min)

Sigh Doubled tidal volume possible every 100 breaths

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003