

**Product Name :**  
Traverse Thread Counter Machines

**Product Code :**  
NLAB-TESTINGLABM65001



**Description :**

Traverse Thread Counter Machines

**Technical Specification :**

Specification:

Fabric Can be Checked: 1/4inch, 1/2 Inch, 1 Inch, 20 mm, 50mm (Calibrated Slots)

Digital Counter: 4-Digit Re-settable Digital

Counter (Registers each increment with a Beep) – 0-9999, L/C: 1

Microscope Magnification: 40x (10x eyepiece and 4x Objective Lens)

Light Intensity Control: Variable light intensity

A fine equipment to check the read/pick of any type of fabric conveniently and accurately

40x magnification with in-built reference line

Calibration Certificate

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003