

**Product Name :**  
Tray Dryers

**Product Code :**  
NLABMEDICAL-M3170004

**Description :**

Tray Dryers India, Tray Dryers manufacturer India, suppliers, Medical Tray Dryers India, Hospital Tray Dryers India

**Technical Specification :**

Gap filled with glass wool for temp. Insulation.

Inner and outer walls made out of thick mild steel sheet on stout mild steel angle with motorized air circulation.

Dial type thermometer has been provided.

Temperature range ambient to 150° C.

Optional:

Digital temp. Indicator-cum-controller.

Digital temp. Indicator-cum-controller with computer data logger

Technical Specifications:

Suitable tray size:

32 x 16 x 1¼"

Chamber sizes in inches	Rating	No. Of trays
33 x 33 x 12	3.0 kw	6
33 x 33 x 17	5.0 kw	12
33 x 33 x 34	7.5 kw	24
68 x 39 x 36	12.0 kw	48
72 x 77 x 36	18.0 kw	96

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003