

Product Name :
Tri Power LED Stereo Microscope

Product Code :
NLAB-EDUSCHOOLEQ173099

Description :

Tri Power LED Stereo Microscope

Technical Specification :

Wide field eyepieces with diopter adjustment on one ocular tube.

Includes a 95 mm black and white reversible stage plate and a frosted glass stage plate.

Rack and pinion focusing with slip clutch and tension adjustment.

Paired objectives: parfocaled, parcentered, and achromatic.

Viewing head has advanced 45° inclined design, and is post mounted to permit 360° rotation.

Adjustable pole mounted style to allow various working distances.

Rugged all metal construction.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003