


<b>Product Name :</b> Tri Power LED Stereo Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173099
--	---



**Description :**

Tri Power LED Stereo Microscope

**Technical Specification :**

Wide field eyepieces with diopter adjustment on one ocular tube.

Includes a 95 mm black and white reversible stage plate and a frosted glass stage plate.

Rack and pinion focusing with slip clutch and tension adjustment.

Paired objectives: parfocaled, parcentered, and achromatic.

Viewing head has advanced 45° inclined design, and is post mounted to permit 360° rotation.

Adjustable pole mounted style to allow various working distances.

Rugged all metal construction.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003