

**Product Name :**  
Triac Characteristics Trainer Analog

**Product Code :**  
NLAB-ELECTRICALAB210017



**Description :**

Triac Characteristics Trainer Analog

**Technical Specification :**

Analog Meters:

Ammeter 500uA DC.

Ammeter 50mA DC.

Voltmeter 3V/30V DC.

Power supplies:

DC Supply IC Regulated 0-5V DC, 150mA.

DC Supply IC Regulated 0-30V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

Triac

Voltage Control through Potentiometer.

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003