

<b>Product Name :</b> Triangular File	<b>Product Code :</b> SWEL-IN-0021NLAB- MACHINELABWK34026
------------------------------------------	-----------------------------------------------------------------



<b>Description :</b>
Triangular File
<b>Technical Specification :</b>
Types: Smooth and Bastard Material: Carbon chrome steel hardened and tempered with rust preventative finish Length: 150mm (6") 200mm (8") 250mm (10")  Standard Welding Equipment Laboratory Set Triangular File Manufacturer,Supplier and Triangular File Exporter in India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003