

Product Name :
Triethylamine AR

Product Code :
NLAB-CHEMICALSIND15001

Description :

Triethylamine AR

Technical Specification :

Triethylamine AR

Molecular Formula: (C₂H₅)₃N

Grade: Analytical Reagents (AR)

Assay: (GC) 99.5%

Capacity: 500ml,2.5Lit,25Lit

Properties: Strong ammonical odour. Toxic by ingestion and inhalation, strong irritant to tissue.

Description: A clear liquid, not more than 20 Hazen units in colour, with an ammonical odour.

Solubility: Miscible with alcohol.

Minimum assay (By GC): 99.5%

Wt. per ml at 20°C: 0.726 - 0.728 g

Refractive Index: 1.400 - 1.402

Water: 0.1%

Lead (Pb): 0.00001%

Non volatile matter: 0.005%

Copper (Cu): 0.000005%

Iron (Fe): 0.00005%

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003