

Product Name : Trinocular Research Microscope	Product Code : NLAB-EDUSCHOOLEQ173066
Description :	
Trinocular Research Microscope	
Technical Specification :	
Specification:- Condenser: Abbe type N A 1.25 with IRIS Diaphram & Filter Eyepiece: WF 10x, 16x, with Eye Guard Objectives: Plan Achromatic 4x, 10x, 40x, 100x Oil Working Stage: Double Layer mechanical stage 140 x 140 mm Focussing Range: Co-axial coarse and fine focussing system Moving Range: 75x50mm. Illumination System: Built in Kohler illumination, Halgon Lamp 6V/20W Viewing Head: Compensation free Trinocular Head	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003