

Product Name :
Triple Beam Balance

Product Code :
NLAB-LBEMALAB230037



Description :

Triple Beam Balance

Technical Specification :

Depending on model, comes with: stainless-steel pan, scoop, 6 in. dia. plate or plastic animal container

Capacity is 2610g, when used with Weight Set (S40037)

Precision ground pivots and patented self-aligning agate bearings to ensure reproducibility

Easy zero adjustment remains firmly in place, while deep-notched, tiered beams with diecast, center-reading poise ensure correct positioning

High-speed throughput allows for more weighings in less time

Beams are easy to read because of their high-contrast markings, while an end-reading device eliminates parallax

Nonmagnetic stainless-steel bearing covers and rust-resistant weighing platforms

Heavy zinc base for exceptional stability, a tamper-resistant poise removable from only one spot on the beam, and attachment weights with recessed pivots to keep them from slipping or falling

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003