

Product Name : Triple Folding Pocket Magnifier	Product Code : NLAB-EDUSCHOOLEQ173086
Description :	
Triple Folding Pocket Magnifier	
Technical Specification :	
A triple-lens magnifier with +5 times magnification for each. This magnifier conveniently folds up for easy transport. Its 2 lenses combine for an impressive magnification to aid in reading or inspecting specimens. Plastic ring approx. 1.25" (30 mm) Dia. Black plastic frame. Magnification 3 lenses: 15X; 2 lenses: 10X; single lens: 5X. Overall length about 1.5" L when closed, 2.5" (65 mm) L when opened.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003