

Product Name :
Triple Spur Gear Chain Pulley Block

Product Code :
NLAB-MECHANICALAB240053



Description :

Triple Spur Gear Chain Pulley Block

Technical Specification :

Available capacity of 500 kg to 10 tons.

Triple spur gear chain pulley block with "ISI" mark & government approved chain pulley block.

"ISI" marked manufacturing ISO 9001:2000 certified company.

Tested as part I0 procedure for 50% over load.

Anti-corrosive powder coated zinc plated parts.

"INDEF" chain pulley blocks made as Per IS 6216 8 8610.

Available in flame proof type.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com