

<b>Product Name :</b> Tripod Magnifier	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173093
<b>Description :</b>	
Tripod Magnifier	
<b>Technical Specification :</b>	
Perfect for viewing rocks, minerals, insects, plants, dissected specimens, etc.	
This is a quality 25mm optical glass 10X lens mounted in a sturdy plastic tripod for steady focusing.	
Tripod size is approximately 3.5 cm tall and 5 cm across top.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003