

Product Name :
Tubing Rubber

Product Code :
NLAB-SCIENTIFICEQ31003



Description :

Tubing Rubber

Technical Specification :

- Medium pressure red colored, extra soft rubber tubing having more wall thickness making it capable of withstanding low to moderately high pressures.
- resistant to acids, alkalis.
- Supplied in Pack of 10m.

Sizes (bore x wall thickness) available are:-

- 3x3.0
- 4.5x3.0
- 6x3.0

- 8x3.0
- 9x3.0
- 10x3.0
- 12.5x4.5

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003