

Product Name :
TVET Automatic Level

Product Code :
NLAB-EDUCATIONTVET40008



Description :

TVET Automatic Level

Technical Specification :

Magnification 32x

Standard Deviation for 1Km Double-Run Leveling : $\pm 1.0\text{mm}$

Objective Aperture : 40mm

Image of Telescope : Erect

Shortest Focusing Distance : 0.65m

Field of View : $1^\circ 20'$

Setting Accuracy : $\pm 0.3''$

Setting Time : 2"

Multiplication Constant : 100

Sensitivity of Circular Level : $8' / 2\text{mm}$

Standard Accessories:

Surveying Tripod ,Adjusting Pin, Desiccating Agent, Hex Key, Plastic Carrying Case, Plumb Bob, User Manual

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003