

Product Name :
TVET Cleveland Open Cup Flash Point Tester

Product Code :
NLAB-EDUCATIONTVET39023



Description :

TVET Cleveland Open Cup Flash Point Tester

Technical Specification :

Power supply: AC 220V \pm 10%, 50 Hz;

Heating device: It is a quartz tube heater, so there is no flame appears.

Explosion protective; its power can be continuously adjusted in the range of 0600W;

Flame application device: automatically scan;

Thermometer: -6~400; Each scale division is 2 ; axles per GB/T514;

Ignition devices:

Ignition source: coal gas (or any other combustible gas);

Spout aperture: 0.60.8mm;

Ambient temperature: -10 ~50;

Relative humidity: 85%;

Total power consumption: not more than 650 W;

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003