

Product Name :
TVET Common Rail Pressure Tester And Simulator

Product Code :
NLAB-EDUCATIONTVET37001



Description :

TVET Common Rail Pressure Tester And Simulator

Technical Specification :

Screen: 84x48 LCD display

Input voltage: 8~32V

Operating current: Typical value <50mA

Operating temperature: 0-50C

Storage temperature: -20C ~70C

Relative humidity: 60%

Product dimensions: 127*72*36

Host machine weight: <250g

Common rail pressure measurement scope: 0-200MPa

Stimulation scope of common rail pressure simulator: 0-200MPa

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003