

**Product Name :**  
TVET Crude Petroleum Water Content Tester

**Product Code :**  
NLAB-EDUCATIONTVET39001



**Description :**

TVET Crude Petroleum Water Content Tester

**Technical Specification :**

Petroleum Products Water Content Tester is designed totally meeting the requirements of GB/T260.

The gripper is designed reasonably. The installation and dismantlement are convenient.

The heating power is continuously adjustable. The voltmeter displays the power institutionally. The temperature control

is reasonable.

It can do parallel test for two samples. The working efficiency is high.

Specifications:

Petroleum Products Water Content Tester (Double Units)

Power supply AC(220±10%)V, 50Hz

Heating power of electric furnace 1000W\*2

Heating control Can be continuously adjusted by a silicon knob

Ambient temperature 35°C

Relative humidity 85%

Maximum power consumption 2200W

Distillation flask: approx. 1000 ml

Receivers: approx. 5ml, graduation is 0.05 ml?

Condenser: 400mm±5mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003