

Product Name :
TVET Distillation Tester

Product Code :
NLAB-EDUCATIONTVET39021



Description :

TVET Distillation Tester

Technical Specification :

Power supply: AC(220±10%)V, 50Hz

Heating power: 1000W,continuous adjustment

Receiving cylinder: 100ml,scale division 1ml

Distillation flask: 125mL. It can meet the requirements GB/T 6536 and ASTM D86

Thermometer: (-2~300)°C and (-2~400)°C. Division value 1°C

Flask support board: SiC,diameter for each hole is 32mm,38mm and 50mm

Ambient temperature: Room temp.~+35°C

Relative humidity: 85%

Maximum power consumption: 1100W

Dimension: 460mm×400mm×500mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003