

Product Name :
TVET Distillation Tower

Product Code :
NLAB-EDUCATIONTVET39026



Description :

TVET Distillation Tower

Technical Specification :

Material stainless steel 304/S316L

Distillation column type packed tower

Theoretical plate 10 ~ 20

Temperature Room temp~200/300oC(depend on the oil heater)

Column size DN50*800-1600mm

Reboiler 20L~30L

Receiver size&quantity 10L*2

Condenser HTA 0.86 (m2)

Vacuity MIN:0.5mbar(depend on the nature of material and vacuum pump&vacuum pump not included)

Floor area length*width*height=1000*800*2500-3300 (mm)

Power 1.1KW(without pump)

Voltage 220V-60HZ

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003