

Product Name :
TVET Drying Oven

Product Code :
NLAB-EDUCATIONTVET40005



Description :

TVET Drying Oven

Technical Specification :

External material: cold-rolled steel with anti-bacteria powder coating.

Stainless steel inner chamber, round angle structure, adjustable shelves.

Double-layer glass observation window, adjustable airtight buckle lock.

PID microprocessor intelligent temperature controller.

Over-temperature protection.

Capacity : 140L

temp. Range : RT+10~250/300?

temp. Precision : 0.1?

temp. Fluctuation : ± 0.5 ?

Timing Range : 1~9999min

Display: LCD

Consumption : ~ 2000 Watt

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003