

Product Name :
TVET Electric Vibrating Sieve Shaker

Product Code :
NLAB-EDUCATIONTVET40013



Description :

TVET Electric Vibrating Sieve Shaker

Technical Specification :

Standard Electric Sieve Shaker for him. 200mm test sieve

Shaking frequency: approx. 1400 times/min.

Shaking method: Up and down

Shaking amplitude: 1.5 mm

Time setting: 0 to 60 min. Selectable

Motor power: 25W

Including Sieve size: Dia. 200 x 50 mm (9pieces)

Sieve opening: 20 to 0.75 mm

Power supply: 220V 50Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003