

**Product Name :**  
TVET Engler Viscometer

**Product Code :**  
NLAB-EDUCATIONTVET39006



**Description :**

TVET Engler Viscometer

**Technical Specification :**

Standard water value:  $51 \pm 1$  s;

Temperature range: 0~100° C;

Temperature controlling accuracy:  $\pm 0.2^\circ$  C;

Engler viscosity thermometer: As per GB/T514;

Graduated flask:  $200 \pm 0.2$ ml;

Inner container: Made by stainless steel;

Accuracy of timer: 0.01s

Timer: Max is 9999s;

Timing mode: timing by a LED;

Heating power: approx. 700W;

Total power consumption: not more than 800 W.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003