

Product Name :
TVET Engler Viscometer

Product Code :
NLAB-EDUCATIONTVET39006



Description :

TVET Engler Viscometer

Technical Specification :

Standard water value: 51 ± 1 s;

Temperature range: $0 \sim 100^{\circ}$ C;

Temperature controlling accuracy: $\pm 0.2^{\circ}$ C;

Engler viscosity thermometer: As per GB/T514;

Graduated flask: 200 ± 0.2 ml;

Inner container: Made by stainless steel;

Accuracy of timer: 0.01s

Timer: Max is 9999s;

Timing mode: timing by a LED;

Heating power: approx. 700W;

Total power consumption: not more than 800 W.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003