

**Product Name :**  
TVET High Intensity Focused Ultrasound

**Product Code :**  
NLAB-EDUCATIONTVET18010



**Description :**

TVET High Intensity Focused Ultrasound

**Technical Specification :**

Superficial dell dermis and SMAS (Superficial Muscular Aponeurotic System).

Fast operation. Full face and neck treatment only takes 20-30 minutes

Completely non-invasive you can put on make-up 15 minutes after treatment.

Instant result will be checked from one time good result will last 2-5 years

Non-invasive painless soft and no side effects

High quality low consumables

More accurate procedure due to Exact Focusing Technology

Power Chip Imported from America--Stable power and Constant depth.

Highest quality Propeller from Korean-- Uniform Point spacing

Equipped 5 heads precisely acting on different depths of the skin

Friendly operating software, adjust energy is OK.

#### TECHNICAL SPECIFICATIONS

Irradiation Type High Intensity Focused Ultrasound

Frequency 4MHz, 7MHz,

Max Energy 2.0J/m<sup>2</sup>

LCD Screen 15 Inch True Color Touch LCD ADSS System

Standard Cartridges 4MHz: 4.5mm, 7MHz: 3.0mm---10,000 shots

Operational Cartridges 2MHz: 13mm(Body/Reduce Fat)

3MHz: 7mm(Body/Reduce Fat)

10MHz: 1.5mm(Eye area)

Face Length: 5~25mm(1mm steps)

Body: 5~25mm(7.0mm, 13mm)

Spot Size Spot<3mm

Dimensions of Machine 45\*30\*15cm

NW/GW 20/22KG

Voltage AC110-220V 50/60Hz

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003