

**Product Name :**  
TVET Hydrafacial

**Product Code :**  
NLAB-EDUCATIONTVET22001



**Description :**

TVET Hydrafacial

**Technical Specification :**

Ultrasound head: promote microcirculation, skin rejuvenating

Multi-polar RF: compact shaping, face lifting

Bio microcurrent: let jaw line tightened, Bio raise the overall tightening of facial skin

Hydrodermabrasion: facial cleansing, peeling, pore cleansing

cold hammer : cold down skin after treatment

oxygen spray gun : skin rejuvenation

Parameters:

Vacuum range 80Kpa

Pressure:600mmHg,Handles:5pcs

Radio frequency 1Mhz, Bipolar

Lon lifting 500Hz (Digital Ion Lifting)

Ultrasound 1Mhz / 2W/cm<sup>2</sup>

User interface 8 inch Color Touch LCD

Input voltage 100~120V, 50/60HZ / 220~240V, 50/60HZ

Noise level 45Db

With 2 ultrasound working mode

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003