

Product Name :
TVET Incubator

Product Code :
NLAB-EDUCATIONTVET41007



Description :

TVET Incubator

Technical Specification :

Intelligent temperature controller, based on microprocessor with PID parameters controlling function, and LCD display make temperature accurate and reliable.

The double door makes the temperature fluctuation small. The toughened-glass inner door is convenient to observe chamber conditions without affecting chamber temperature.

The working chamber adopts 304 stainless steel and the round angle structure , easy to clean.

Capacity : 50 Liters

temp. Range : RT+5~65?

Time Setting : 1~9999mins

Alarm : Over-temperature alarm; audible and visual alarm.

Shelves : 2 pcs/adjustable

Power Consumption : ~250 Watt

Power Supply : AC110/220V±10%, 50/60Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003