

Product Name :
TVET Infrared Sauna Heater Room

Product Code :
NLAB-EDUCATIONTVET22003



Description :

TVET Infrared Sauna Heater Room

Technical Specification :

Dimensions : Approx. 1800x1500x2100 mm

GW : 178kgs

Capacity : 4 people

Material : Canadian Hemlock or red cedar

Wall Construction: Double-wall construction

Safety Device: Double Insulation /Infrared

Power : 2700W

Voltage : 110V/220V/380V

Heater : Far infrared carbon panel

Frequency : 50HZ/60HZ

Door Glass : 8mm tempered glass

Control panel : Inside and outside control panel (digital with USB)

Handle : Wooden handle and Luxury crystal handle for option

Accessories : Oxygen ionizer, USB radio, Speakers, Lamps

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003