

Product Name :
TVET LCD Digital Biological Microscope

Product Code :
NLAB-EDUCATIONTVET39009



Description :

TVET LCD Digital Biological Microscope

Technical Specification :

Viewing Head: Compensation Free Binocular Head, Inclined at 30°, Interpupillary Distance 48-75mm

Camera

System Valid Pixel

2048x1036(3.0M Pixel)

Output Mode USB2.0

Operation System WINDOWS7/8 2000 / XP/ VISTA

Software ScopelImage 9.0

Range of Viewing Field 90%

Eyepiece: WF10x/ 18

Objective:

Infinite planAchromatic Objective 4x,10x,40x,100x

Nosepiece:

Quadruple Nosepiece

Stage:

Double Layers Mechanical Stage 140x140mm Moving range 75x40mm

Condenser:

Abbe NA1.25 sliding-in centerable condenser

Focusing: Coaxial Coarse & Fine Adjustment System, Range 24mm, Fine Division 0.004mm

Coaxial Coarse & Fine Adjustment System, Range 20mm,

Illumination:

220V or 110V/6V/20W Halogen Lamp, Brightness Adjustable.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003