

Product Name :
TVET Light Therapy

Product Code :
NLAB-EDUCATIONTVET18003



Description :

TVET Light Therapy

Technical Specification :

Adopts LED light therapy to accelerate the skin circulation for anti-aging, no side effect, healthy and safe.

6 colors light, which has different wave length for light treatment of skin problem.

It can effectively reduce wrinkles and sagging, shrink pores, tighten and firm skin, make skin tender, elastic and young.

Improve dark complexion caused by bad metabolism or poor blood circulation.

Repair and nurse the damaged skin, easy to operate, comfortable and painless.

A good skin care tool for solving skin problems, such as yellow and dark skin, acne, large pores, loose skin, etc.

Function

Red(640nm):Deep repair and increase skin elasticity.

Blue(415nm):Remove acne and inhibit bacteria.

Yellow(590nm):Improves wrinkle and enhances immunity.

Green(532nm):Lighten spots and calming skin.

Purple: Improves fine lines

Cyan:Deep metabolism and control oil.

Specification

Item Type:PDT

Facial Mask

Color:White

Material:ABS

Input:100-240V 50/60Hz

Output:24V 5A

Weight:1835g

Size:45*25*31cm/17.72*9.84*12.2in

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003