

Product Name :
TVET Liquid Nitrogen Tanks

Product Code :
NLAB-EDUCATIONTVET22009



Description :

TVET Liquid Nitrogen Tanks

Technical Specification :

Features of large capacity liquid nitrogen container:

Made by high strength aluminum alloy.

Can use draw the square canister, easy to identify and access sample independently.

Can storage human, animal organs, install for large size workpiece.

High vacuum multi-layer heat insulation design, no less than 5 years vacuum qualify

Volume 60L

Neck diameter 210mm

External diameter 560mm

Height 900mm

Empty weight 32kg

Liquid nitrogen storage time 63-74days

Liquid nitrogen evaporation 0.851-1.0L/day

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003