

**Product Name :**  
TVET Liquid Petroleum Products Hydrocarbon Tester

**Product Code :**  
NLAB-EDUCATIONTVET39003



**Description :**

TVET Liquid Petroleum Products Hydrocarbon Tester

**Technical Specification :**

Power supply: AC220V±5%,50Hz

Pressure regulating range of reducing valve:(0~400)kPa

Electric agitator: independently controlled for each way

Ultraviolet light source pipe: 1220mm in length,wavelength is 365mm±5nm

Illuminating lamp: 1220mm in length, power is 40W

Air supply: Nitrogen cylinder (or air compressor, compressed air bottle)

Ambient temperature: (5~35)°C

Relative humidity: 85%

Dimension: 350mm×400mm×1770mm

Maximum power consumption: 600W

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003