

Product Name :
TVET Magnetic Hotplate Stirrer

Product Code :
NLAB-EDUCATIONTVET41005



Description :

TVET Magnetic Hotplate Stirrer

Technical Specification :

Work plate material: Glass,ceramic

Maximum Power [W]: 1030

Stirring positions: min.1

Stirring quantity [H2O]: 10 L

Max. magnetic bar length [mm]: 80

Speed range [rpm]: 0 ~ 1500

Heating output [W]: 1000

Heating temperature range [°C]: RT~550°C,increment 5

Safety temperature [°C]: 580

Power Supply : AC110/220V±10%, 50/60Hz

Work Plate Size : 184*184mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003