

**Product Name :**  
TVET Maintenance Tool Bag

**Product Code :**  
NLAB-EDUCATIONTVET20003



**Description :**

TVET Maintenance Tool Bag

**Technical Specification :**

5 long nose pliers.

4 diagonal pliers.

6 lineman pliers.

Desoldering braids.

Solder cores.

Desoldering pump.

Snap-off blade cutters.

Stainless tweezers.

Needle file round and flat 140mm.

Normal grip screwdriver:

Slotted : 3x5mm, 4x6mm.

Phillips : #1x3x3.2mm, #1x3x5mm.

Nuts : 3/16x3.5.

High-speed soldering iron 20W-130W switchable.

7 hhok-scraper.

7 fork spikes.

7 brush-knife.

Alignment screw driver.

Solder paste.

Carrying zip bag.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003