

**Product Name :**  
TVET Manual Lathe Machine

**Product Code :**  
NLAB-EDUCATIONTVET20004



**Description :**

TVET Manual Lathe Machine

**Technical Specification :**

Swing over bed : 355 mm

Swing over carriage : 220 mm

Swing over gap : 500 mm

Width of bed-way : 186 mm

Distance between centers : 1000 mm

Taper of spindle : MT5

Spindle diameter: 38 mm

Range of speed : 70~2000 rpm

Step of speed : 8

Head : D1-4

Metric thread : 26kinds (0.4~7mm)

Inch thread : 34kinds(4~56T.PI)

Longitudinal feeds (mm/r) : 16 kinds(0.35~5M.P)

Cross feeds (mm/r) : 36 kinds(6~104D.P)

Standard accessories:

Three jaw chuck, Four jaw chuck, Steady rest, Follow rest, Working light, Coolant system, Foot brake, Tools and tool box & Instruction Manual

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003