

Product Name :
TVET Manual Strain Direct Shear Apparatus

Product Code :
NLAB-EDUCATIONTVET40009



Description :

TVET Manual Strain Direct Shear Apparatus

Technical Specification :

Area of the soil sample: 30cm²

Sample thickness: minimum 2 cm

Lever ratio: 1:12

Vertical loading: 50, 100, 200, 300, 400 (Kpa)

Max. horizontal loading: 1.2KN

Relative weight of weights: approx. 1.275, 2.55, 5.1, 7.65, 10.2 Kpa, loading hanger axles per grade 1

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

