


**Product Name :**  
TVET Marsh Funnel

**Product Code :**  
NLAB-EDUCATIONTVET39015



**Description :**

TVET Marsh Funnel

**Technical Specification :**

Mesh size 1.6mm(12 mesh)

The volume below the crunchy 1500ml±15ml

Precision When 1500ml standard distilled water was injected into the funnel, the outflow time of 946ml standard distilled water is 26±0.5s

Overall dimensions 570×157×362mm

Packing size 385×205×230mm

## Naugalabequipments

**Website:** www.naugalabequipments.com, **Email:** sales@naugalabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003