

Product Name :
TVET Marsh Funnel

Product Code :
NLAB-EDUCATIONTVET39017



Description :

TVET Marsh Funnel

Technical Specification :

Mesh filter diameter 1.6mm(12 mesh)

Funnel Capacity 1500±15ml

Accuracy When to funnel the standard of 1500ml of distilled water,

out of 946ml standard distilled water in time for the 26±0.5s

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003