

Product Name :
TVET Microscope

Product Code :
NLAB-EDUCATIONTVET39007



Description :

TVET Microscope

Technical Specification :

Viewing Head :Sliding Trinocular Head Inclined at 45° -Nosepiece

Ball Beating Quadruple Nosepiece

Eye-piece: WF10X/18mm

Objective Achromatic :4X,10X,40X(S)100X(S,O)

Stage Double Layer Mechanical Stage 135*140mm

Condenser ABBE NA1. 25 condenser with Iris Diaphragm &filter,rack&pinion adjustable

Illumination Build-in Illumination,Halogen lamp6V/20W,Power supply 90-230V

PowerSupply AC110/220V±10%,50/60Hz

Optional Accessory Camera, Monitor

Package Cartoon with foam

Package Size 250*360*460mm

Gross Weight 8kg

Image Sensor 1/3 SONY CCD Video Output BNC 1.0VP-P.75

-Signal System PAL/NTSC Backlight Compensation OFF/ON

ScanningFrequency Horizontal:15.62KHz

Vertical:50Hz Electronic Shuttle AUTO IRIS EE

Horizontal Resolution 480TV/520TV Camera Interface CS

Min Illumination F1.2 0.5LUX

0.01 LUX Auto Balance

Gain Control Auto Power DC 12V

48S/N Ratio 50dB AGC

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003