

Product Name :
TVET Milking Machine

Product Code :
NLAB-EDUCATIONTVET41010



Description :

TVET Milking Machine

Technical Specification :

Milking Trolley

Pulsation mode: pneumatic pulsation

Pulsation beat: single beat or binaural beats

Pulsation times: 64~75 times/min

Operation vacuum degree: approx. 50Kpa

Production capacity: approx. 20~24 heads/h

Vacuum Pump

Air Flow: approx. 12m³/h

Power: approx. 1.1KW,220V

Rotary Speed: approx. 1400r/min

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003