

Product Name :
TVET Milkshake Machine

Product Code :
NLAB-EDUCATIONTVET22007



Description :

TVET Milkshake Machine

Technical Specification :

Zinc Motor Housing and base with stainless steel brushed appearance;

Stainless Steel Cup with approx. 0.9 liter capacity;

Stepless speed regulation switch or high/low switch;

Voltage/Frequency: 220(V)/50(Hz)

Power: approx. 0.3 (kW)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003