

**Product Name :**  
TVET Milkshake Machine

**Product Code :**  
NLAB-EDUCATIONTVET22007



**Description :**

TVET Milkshake Machine

**Technical Specification :**

Zinc Motor Housing and base with stainless steel brushed appearance;

Stainless Steel Cup with approx. 0.9 liter capacity;

Stepless speed regulation switch or high/low switch;

Voltage/Frequency: 220(V)/50(Hz)

Power: approx. 0.3 (kW)

## Naugalabequipments

**Website:** www.naugalabequipments.com, **Email:** sales@naugalabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003