

**Product Name :**  
TVET Model Sedimentation Tank

**Product Code :**  
NLAB-EDUCATIONTVET28002



**Description :**

TVET Model Sedimentation Tank

**Technical Specification :**

This equipment has been designed to study the hydraulic behaviour and the settling efficiency of a sedimentation tank.

Framework of AISI 304 stainless steel with castors

Tank of clear methacrylate with capacity of 80 litres and section for mixing tap water and suspension

Suspension feed tank with capacity of 120 litres and recycle

Centrifugal pump of AISI 304 stainless steel

Variable area flowmeter for tap water

Variable area flowmeter for suspension

Power supply: 230 Vca 50 Hz single-phase - 0.5 KVA

(other voltage and frequency on demand)

Dimensions: 1700 × 800 × 1500 (h) mm

Weight: 120 Kg

Supplied theoretical-experimental manual

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003