

Product Name :
TVET Multishape Machine

Product Code :
NLAB-EDUCATIONTVET18001



Description :

TVET Multishape Machine

Technical Specification :

- Wrinkle removal
- skin lift
- Skin tightening
- Skin surface smooth
- Eyelid treatment area
- Skin surface smooth
- Body shaping
- Body slimming
- Cellulite reduction
- Massage

Body circumference reduction

Technical Data :

Input Power: 220V; 50Hz/60Hz

Cavitation Frequency: 1MHz, 50W(BLT)

Vacuum: 680mmHg, 1~15Hz Pulse(Continuous/ Pulse/ Super Pulse)

RF Frequency: Multi-polar RF

Power 50W with light therapy

Power 50W with pulse

Display 10.2" Color Touch LCD

Light Therapy Red LED (630nm), Blue LED(430nm)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003