

**Product Name :**  
TVET Pellet Mixer

**Product Code :**  
NLAB-EDUCATIONTVET41011



**Description :**

TVET Pellet Mixer

**Technical Specification :**

Capacity 300 Kg/hour, 25 kg/process

Productivity Time : 10 ~ 15 Minutes

Power : Diesel 16 HP

Dimensions: 110 x 90 x 135 cm

Material : Stainless Steel Blade, Steel Frame/Housing

Frame : Steel

Transmission : Pulley

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com