

Product Name :
TVET PIC Trainer

Product Code :
NLAB-EDUCATIONTVET26002



Description :

TVET PIC Trainer

Technical Specification :

Training box main parameters:

Overall Size: 420 (L) x 300 (W) x 120mm (H); 3.5KG

The power supply: single-phase power supply, AC 220 V, 50/60 Hz, power <50VA.

With grounding protection, over-current protection, power supply indication, main power indication

Working environment: temperature: -10?~+40?

Signal Source Technology Parameter:

DC 5V voltage source, capacity: 1 A

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com