

Product Name :
TVET Polish Polisher Machine

Product Code :
NLAB-EDUCATIONTVET19007



Description :

TVET Polish Polisher Machine

Technical Specification :

Power: approx. 1100W

Voltage: 220V

Base plate diameter: approx. 14"

Speed: approx. 140rpm/min

Main cable length: approx. 8 m

Multi functions: floor cleaning, polishing, waxing

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com