

Product Name :
TVET Polished Rod

Product Code :
NLAB-EDUCATIONTVET40011



Description :

TVET Polished Rod

Technical Specification :

Type: polished rod

External diameter

(+0.127,-0.254)1 1/8(28.6)

Thread pin nominal value (mm): 15/16(23.8)PR,1 1/16(27)PR

Matched sucker rod (mm): 5/8(15.9),3/4(19.1)

Mechanical properties: carbon steel?620MPa? tensile strength?827MPa

steel?655MPa?tensile strength?1100MPa

Coating properties: coating hardness 480HV200

coating thickness 0.254~0.508mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003