

Product Name :
TVET Polishing Machine

Product Code :
NLAB-EDUCATIONTVET40002



Description :

TVET Polishing Machine

Technical Specification :

Automatic grinding input: 0.04- 0.12mm

Grinding spec. Specimens: 50-150mm cubes dia. 50-dia.120mm cylinders

Degree of finish (by standards of rough machining): Ra 3.2

The precision of standard flatness: ± 0.05 mm (within 50 range)

Diameter of the grinding wheel: 200mm

Cooling devices: receive water cooling

Voltage: 380V, three-phase

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com